

DIMENSIONS OF PDL-COATED GLASS-BOTTOM VIEWPLATE-96 F (PN 6005531, 6005530)

Material

Polystyrene plate frame with black pigment. Transparent glass bottom. Material is Tissue Culture treated and poly-D-lysine coated to facilitate cell adherence. Product is sterilized by gamma irradiation.

Dimensions

Length of plate	127.75	mm
Width of plate	85.34	mm
Height of plate without lid	11.60	mm
Plate bottom thickness	0.175	mm
Distance of bottom to ground	0.50	mm
A1 offset to top	11.30	mm
A1 offset to side	14.30	mm
Well depth	11.35	mm
Working volume	50 – 200	µL
Max. well volume	300	µL

