


DIMENSION OF UNIFILTER-96 MICROPLATES (PN 6005174, 6005177)

Material

White Barex microplate with 96-shallow wells (150 µm) with glass fiber bottoms supported by plastic grid.

Dimensions

Length of plate	127.80	mm
Width of plate	85.60	mm
Height of plate without lid	14.60	mm
Tolerances of dimensions	± 0.13	mm
Filter thickness	GF/C 0.68	mm
	GF/B 0.26	mm
A1 offset to top	11.30	mm
A1 offset to side	14.30	mm
Well depth	3.20	mm
Working volume	50 – 100	µL
Max. well volume	150	µL



Product distributed by:
PerkinElmer LAS B.V.
The Netherlands