


DIMENSIONS OF CELLCARRIER-96 MICROPLATES (PN 6005550, 6005920)



Total volume per well: 392 μ L

Cat#	Description	A	B	C	D	D'	E	F	G	H	I
6005550	CellCarrier-96	127.76 ± 0.2	85.48 \pm 0.2	14.4 ± 0.1	\varnothing 6.96 0.05 -0.15	\varnothing 6.58 0.05 -0.15	2.5 \pm 0.05	10.9 ± 0.1	14.38 ± 0.2	11.24 ± 0.2	9 \pm 0.05
6005920	Collagen Coated CellCarrier-96	127.76 ± 0.2	85.48 \pm 0.2	14.4 ± 0.1	\varnothing 6.96 0.05 -0.15	\varnothing 6.58 0.05 -0.15	2.5 \pm 0.05	10.9 ± 0.1	14.38 ± 0.2	11.24 ± 0.2	9 \pm 0.05

All dimensions in mm