


DIMENSION OF UNIFILTER-96 MICROPLATES (PN 6005174, 6005177)

Material

White Barex microplate with 96-shallow wells (150 µm) with glass fiber bottoms supported by plastic grid.

Dimensions

| | | |
|-----------------------------|-----------|----|
| Length of plate | 127.80 | mm |
| Width of plate | 85.60 | mm |
| Height of plate without lid | 14.60 | mm |
| Tolerances of dimensions | ± 0.13 | mm |
| Filter thickness | GF/C 0.68 | mm |
| | GF/B 0.26 | mm |
| A1 offset to top | 11.30 | mm |
| A1 offset to side | 14.30 | mm |
| Well depth | 3.20 | mm |
| Working volume | 50 – 100 | µL |
| Max. well volume | 150 | µL |



Product distributed by:
PerkinElmer LAS B.V.
The Netherlands