

## **Phosphatase Labeled Antibody to Mouse IgG (H+L) (HAS) [Goat]**

**Catalog #:** NEF824001EA

**Size:** 1.0 ml

**Description:** Affinity purified antibody, isolated from a pool of serum goats immunized with purified Mouse IgG was labeled with phosphatase using a modified glutaraldehyde procedure.

**Form / Storage:** Liquid. Store @ 2 to 8°C

**Contents:** Each vial contains 1.0 ml @ 1.0 mg/ml of affinity purified Goat Anti-Mouse IgG .

**E/P Ratio:** Molar enzyme / antibody protein ration = 2.7 : 1.

**Specificity:** Tested by immunoelectrophoresis, gel diffusion & ELISA techniques as applicable. This product reacts specifically with Mouse IgG and may recognize other immunoglobulin types that have light chains in common with IgG. Cross-reactivity to human serum has been minimized using affinity procedures. No antibody is detectable to non-immunoglobulin serum components.

**Stabilizers:** Stabilization diluents added as a protein stabilizer. No preservatives added.

**Stability:** This product is stable for one year from date of receipt when stored properly.

PerkinElmer, Inc.  
940 Winter Street  
Waltham, MA 02451 USA  
Phone: (800) 762-4000 or  
(+1) 203-925-4602  
[www.perkinelmer.com](http://www.perkinelmer.com)

For a complete listing of our global offices, visit [www.perkinelmer.com/contactus](http://www.perkinelmer.com/contactus)